

# Q S Series - Sharps Containers



**SANI-Safe**  
200ml - : QSsn



**RIGI-Box mini**  
300ml - : QSrigi



**QUICK-Safe 0.5**  
500ml - : QS 0.5



**QUICK-Safe 1.0**  
1.0 Litre - : QS 1.0



**OPTI-Safe**  
700ml, 1.7, 2.1 & 3.0 Litre  
QSopt 0.7, 1.7, 2.1, 3.0



**VARIO**  
1.5 & 2.3 Litre  
QSvo 1.5, 2.3



**EURO-Safe**  
3.0, 4.0, 6.0, & 7.0 Litre  
QSeu 3.0, 4.0, 6.0, 7.0



**STERI-Safe**  
35 & 60 Litre - : QSsi 35, 60

**Medical Disposal Containers**

**Idc medical**

16/13 Gibbens Road, West Gosford  
NSW Australia 2250  
ACN 078 760 504 ABN 84 078 760 504

Phone: 02 4321 0299  
Internet: [www.idcmedical.com](http://www.idcmedical.com)

Fax: 02 4321 0288  
eMail: [sharps@idcmedical.com](mailto:sharps@idcmedical.com)